
RESORBA® Barbed Suture

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- siliconised needles made of 300 steel
- synthetic with a long period of action
- unidirectional spiral angle cut on one side
- high tensile strength
- triangular end stopper

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RESORBA® Barbed Suture is a polydioxanone knotless tissue-closure device, which is comprised of barbed dyed polydioxanone (PDO) suture material, armed with a surgical needle on one end and end stopper at the other end.

The RESORBA® Barbed Suture bars are unidirectional oriented to allow for tissue approximation without the need to tie surgical knots.

RESORBA® Barbed Suture is consisting of polyester, poly(p-dioxanone) and is dyed with D&C violet No. 2. While the formation of barbs in the RESORBA® Barbed Suture reduces the tensile strength relative to non-barbed suture material of the same size, tying of knots in non-barbed suture materials also reduces their effective strength. For this reason, the strength of the RESORBA®

Barbed Suture can be compared to USP knot strength of non-barbed sutures.

After implantation the breaking load is at 55% after 42 days. Absorption by hydrolysis is completed after about 180 – 220 days.

- Colour: violet
- Chemical name: poly (p-dioxanon)
- Thread diameter: USP 4-0 -1 (1.5-4 metric)
- Types of packaging: needle-thread-combinations
- Sterilization method: ethylene oxide

USP strength before barbed suture	USP strength of RESORBA® Barbed Suture	Number of barbs per inch
2	1	32
1	0	33
0	2-0	42
2-0	3-0	51
3-0	4-0	54

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Properties of RESORBA® Barbed Suture

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Needle


300 series steel needle provides high bending stiffness and surface siliconisation improves glide through the tissue.

Barbed Suture

The number of **unidirectional barbs** is adapted to the respective thread strength. The optimized distance between the individual barbs and their special shape ensure secure and firm fixation across all thread sizes.

End stopper

The triangular end stopper eliminates the need for "looping" at the beginning of suturing, making handling quicker and easier. The end stopper is perpendicular to the tissue and thus prevents the suture from running out and fixes the suture in the desired position.

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


Holding the needle

Needles should be held approx. 3/4 away from the tip. Do not clamp where the suture is swaged to the needle to avoid weakening the needle and suture.

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
Straightening the suture

Straighten the suture gently after removal. Do not rub it or pull it quickly. Do not hold the needle and stretch!

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1st letter: Indicates the curvature of the needle

F = 5/8 circle
H = 1/2 circle
D = 3/8 circle
V = 1/4 circle
K = semicurved
A = asymptotic
G = straight

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2nd letter: Gives information on the type of needle and needle tip

R = round bodied needle
S = reverse cutting needle

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3rd and 4th letters: escribe special features of the needle and needle tip

I = conventional cutting
M = PREMIUM cut (partly hand-honed)
N = blunt point
T = trocar point
P = spatulated needle
S = stronger diameter
X = extra strong diameter
F = thin PREMIUM cut "THIN LINE"
(partly hand-honed)
W = flexible needle

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Numbers indicate the straight (extended) length of the needle in millimetres

S (after numbers) = stronger diameter
F (after numbers) = extra thin diameter

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
1. Spatula needle □ = P

1/2 circle, 3/8 circle, 1/4 circle, or straight = HSPM, DSPM, VSPM

- For use in ophthalmic surgery and microsurgery
- Needle cross-section flattened
- PREMIUM cut
- Cuts on lateral edges

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2. Blunt round-bodied needle ○ = RN

1/2 circle, 3/8 circle or semicurved = HRN, DRN, KRN

- For use in parenchymal tissue, the cervix, and muscle traction in the eye
- Needle tip blunt
- No puncturing of vessels or tendons

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
3. Reverse cutting needle ▼ = S

1/2 circle, 3/8 circle, semicurved or straight = HS, DS, KS, GS

- For use in firm tissue such as skin
- Triangular needle cross-section
- Also available with PREMIUM cut needles M

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
5. Round-bodied needle ● = R

5/8 circle, 1/2 circle, 3/8 circle or straight = FR, HR, DR, GR

- For use in soft tissue (subcutaneous) such as muscle, fascia, mucosa
- In order to improve grip by the needle holder, the cross-section of the middle part of the needle is flattened
- Easy insertion

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▽ = RT

1/2 circle, 3/8 circle or straight = HRT, DRT, GRT, ART

- For use in firm tissue, sclerotic vessels, and implants
- Needle tip triangular in section

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Absorbable sutures

approximate the tissues during the healing process. During this time the suture's tensile strength will gradually diminish. Absorbable suture material is metabolised by endogenous proteolytic enzymes or by hydrolysis (in the case of PGA **RESORBA®**, PGA **resoquick™**, PDO **RESORBA™** and **GLYCOLON®**).

Non-absorbable sutures

remain almost unchanged when placed within body tissues and is encapsulated within the wound scar tissue by the organism. The sutures used for skin closure are removed once the scar tissue has become sufficiently firm to

hold the wound edges together (usually after 7 – 14 days).

It has to be distinguished as follows:

Absorption time

The period in which the suture loses 50 % of its knot tensile strength.

Disintegration


The period during which essentially non-absorbable suture break down by degradation into (smaller) pieces and thus loosing its strength (e.g. polyamide).

Dissolution

The period during which the suture completely dissolves within the tissue.

Characteristics of absorption

Different indications also require different tensile strengths and absorption characteristics. These particular features of different sutures can be achieved by the choice of the material and modifying the production process. In addition to the immediate, moderately quick or delayed loss of tensile strength there is also the corresponding duration of absorption. Any given thread material can only fulfil its purpose as long as it has the desired tensile strength.

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*In-vitro trial of suture degradation by measuring the knot tensile strengths of **GLYCOLON®**, **PGA resoquick™**, **PGA RESORBA®** and **PDO RESORBA™**. Suture size: 3-0 USP (2 metric).*

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RESORBA Medical GmbH, Am Flachmoor 16, 90475 Nürnberg, Germany