
RESORBA® Barbed Suture

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- *siliconised needles made of 300 steel*
- *synthetic with a long period of action*
- *unidirectional spiral angle cut on one side*
- *high tensile strength*
- *triangular end stopper*

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RESORBA® Barbed Suture is a polydioxanone knotless tissue-closure device, which is comprised of barbed dyed polydioxanone (PDO) suture material, armed with a surgical needle on one end and end stopper at the other end. The RESORBA® Barbed Suture barbs are unidirectional oriented to allow for tissue approximation without the need to tie surgical knots.

RESORBA® Barbed Suture is consisting of polyester, poly(p-dioxanone) and is dyed with D&C violet No. 2. While the formation of barbs in the RESORBA® Barbed Suture reduces the tensile strength relative to non-barbed suture material of the same size, tying of knots in non-barbed suture materials also reduces their effective strength. For this reason, the strength of the RESORBA®

Barbed Suture can be compared to USP knot strength of non-barbed sutures. After implantation the breaking load is at 55% after 42 days. Absorption by hydrolysis is completed after about 180 – 220 days.

- ➔ Colour: violet
- ➔ Chemical name: poly (p-dioxanone)
- ➔ Thread diameter: USP 4-0 -1 (1.5-4 metric)
- ➔ Types of packaging: needle-thread-combinations
- ➔ Sterilization method: ethylene oxide

USP strength before barbed suture	USP strength of RESORBA® Barbed Suture	Number of barbs per inch
2	1	32
1	0	33
0	2-0	42
2-0	3-0	51
3-0	4-0	54

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Properties of RESORBA® Barbed Suture

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Needle


300 series steel needle provides high bending stiffness and surface siliconisation improves glide through the tissue.

Barbed Suture

*The number of **unidirectional barbs** is adapted to the respective thread strength. The optimized distance between the individual barbs and their special shape ensure secure and firm fixation across all thread sizes.*

End stopper


The triangular end stopper eliminates the need for "looping" at the beginning of suturing, making handling quicker and easier. The end stopper is perpendicular to the tissue and thus prevents the suture from running out and fixes the suture in the desired position.

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
Short lesson in symbols used for medical products

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
Holding the needle

Needles should be held approx. 3/4 away from the tip. Do not clamp where the suture is swaged to the needle to avoid weakening the needle and suture.

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
Straightening the suture


Straighten the suture gently after removal. Do not rub it or pull it quickly. Do not hold the needle and stretch!

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
1. Spatula needle ◻ = P

1/2 circle, 3/8 circle, 1/4 circle, or straight = HSPM, DSPM, VSPM

- For use in ophthalmic surgery and microsurgery
- Needle cross-section flattened
- PREMIUM cut
- Cuts on lateral edges

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
2. Blunt round-bodied needle ○ = RN

1/2 circle, 3/8 circle or semicurved = HRN, DRN, KRN

- For use in parenchymal tissue, the cervix, and muscle traction in the eye
- Needle tip blunt
- No puncturing of vessels or tendons

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
3. Reverse cutting needle ▼ = S

1/2 circle, 3/8 circle, semicurved or straight = HS, DS, KS, GS

- For use in firm tissue such as skin
- Triangular needle cross-section
- Also available with PREMIUM cut needles M

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
5. Round-bodied needle ● = R

5/8 circle, 1/2 circle, 3/8 circle or straight = FR, HR, DR, GR

- For use in soft tissue (subcutaneous) such as muscle, fascia, mucosa
- In order to improve grip by the needle holder, the cross-section of the middle part of the needle is flattened
- Easy insertion

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◐ = RT

1/2 circle, 3/8 circle or straight = HRT, DRT, GRT, ART

- For use in firm tissue, sclerotic vessels, and implants
- Needle tip triangular in section

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Absorbable sutures

*approximate the tissues during the healing process. During this time the suture's tensile strength will gradually diminish. Absorbable suture material is metabolised by endogenous proteolytic enzymes or by hydrolysis (in the case of PGA *resoquick*TM, GLYCOLON[®], PGA *RESORBA*[®], PDO *RESORBA*TM and *RESORBA*[®] Barbed Suture).*

Non-absorbable sutures

remain almost unchanged when placed within body tissues and is encapsulated within the wound scar tissue by the organism. The sutures used for skin closure are removed once the scar tissue

has become sufficiently firm to hold the wound edges together (usually after 7 – 14 days).

It has to be distinguished as follows:

Absorption time

The period in which the suture loses 50 % of its knot tensile strength.

Disintegration


The period during which essentially non-absorbable suture break down by degradation into (smaller) pieces and thus losing its strength (e.g. polyamide).

Dissolution

The period during which the suture completely dissolves within the tissue.

Characteristics of absorption

Different indications also require different tensile strengths and absorption characteristics. These particular features of different sutures can be achieved by the choice of the material and modifying the production process. In addition to the immediate, moderately quick or delayed loss of tensile strength there is also the corresponding duration of absorption. Any given thread material can only fulfil its purpose as long as it has the desired tensile strength.

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*The graph shows a simplified representation of the respective tear force losses over time. Data of GLYCOLON[®] and PGA *RESORBA*[®] come from in-vitro experiments, whereas the data of PDO *RESORBA*TM and PGA *resoquick*TM come from in-vivo experiments. Thus, the curves are not directly comparable but serve as a simplified visual representation.*

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