
RESORBA® Barbed Suture

RESORBIERBARES NAHTMATERIAL

resorba.com

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
- ✓ *siliconised needles made of 300 steel*
- ✓ *synthetic with a long period of action*
- ✓ *unidirectional spiral angle cut on one side*
- ✓ *high tensile strength*
- ✓ *triangular end stopper*

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► RESORBA® Barbed Suture (PDO) is a polydioxanone knotless tissue-closure device, which is comprised of barbed dyed polydioxanone (PDO) suture material, armed with a surgical needle on one end and an end stopper at the other end.

RESORBA® Barbed Suture (PDO) features unidirectionally oriented barbs that facilitate tissue approximation without the need for tying surgical knots.

Although the barbs in the RESORBA® Barbed Suture (PDO) reduce tensile strength compared to non-barbed

suture material of the same size, tying knots in non-barbed sutures also diminishes their effective strength. For this reason, the strength of the RESORBA® Barbed Suture (PDO) can be compared to USP knot strength of non-barbed sutures.

After implantation, the tensile strength remains at 55 % after 42 days. Absorption by hydrolysis is completed after about 180 to 220 days. RESORBA® Barbed Suture (PDO) is available in violet (coloured with the physiologically harmless D+C No. 2 dye).

USP strength before barbed suture	USP strength of RESORBA® Barbed Suture	Number of barbs per inch
2	1	32
1	0	33
0	2-0	42
2-0	3-0	51
3-0	4-0	54

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Properties

- ✓ siliconised needles made of 300 steel
- ✓ synthetic with a long period of action
- ✓ unidirectional spiral angle cut on one side
- ✓ high tensile strength
- ✓ triangular end stopper

- **Colour:** violet
- **Chemical name:** poly(p-dioxanone)
- **Thread diameter:** USP 4-0 - 1 (1.5-4 metric)
- **Types of packaging:** needle-thread-combinations
- **Sterilization method:** ethylene oxide

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Needle

The 300-series steel needle offers exceptional bending stiffness, while its siliconized surface ensures smooth tissue penetration and optimal glide.


Barbed Suture

*The number of **unidirectional barbs** is adapted to the respective thread strength. The optimized distance between the individual barbs and their special shape ensure secure and firm fixation across all thread sizes.*

End stopper

The triangular end stopper eliminates the need for "looping" at the beginning of suturing, making handling quicker and easier.


The end stopper is perpendicular to the tissue and thus prevents the suture from running out and fixes the suture in the desired position.

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


The thread must be stretched gently after it has been removed from the primary packet. Do not pull or rub it

abruptly. Do not grasp the needle and stretch the thread!

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Needles should be held approx. 3/4 away from the tip. Do not clamp where the suture is swaged to the needle to avoid weakening the needle and suture.

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1. letter: Indicates the curvature of the needle

- F = 5/8 circle
- H = 1/2 circle
- D = 3/8 circle
- V = 1/4 circle
- K = semicurved
- J = J-shaped
- A = asymptotic
- T = compound circle
- G = straight

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2. letter: Gives information on the type of needle and needle tip

- R = round bodied needle
- S = reverse cutting needle

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3.+4. letters: Describe special features of the needle and needle tip


- I = conventional cutting
- M = PREMIUM cut (partly hand-honed)
- N = blunt point
- T = trocar point
- P = spatulated needle
- S = stronger diameter
- X = extra strong diameter
- F = thin PREMIUM cut "THIN LINE" (partly hand-honed)
- W = flexible needle

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
Numbers indicate the straight (extended) length of the needle in mm

- S (after number) = stronger diameter
- F (after number) = extra thin diameter

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
Spatula needle ◻ = P

1/2 circle, 3/8 circle, 1/4 circle, or straight
= HSPM, DSPM, VSPM

- ▶ For use in ophthalmic surgery and microsurgery
- ▶ Needle cross-section flattened
- ▶ PREMIUM cut
- ▶ Cuts on lateral edges

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Reverse cutting needle ▼ = S

1/2 circle, 3/8 circle, semicurved or straight
= HS, DS, KS, GS

- ▶ For use in firm tissue such as skin
- ▶ Triangular needle cross-section
- ▶ Also available with PREMIUM cut needles M

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

Round-bodied cutting needle ◐ = RT

1/2 circle, 3/8 circle, J-shaped, asymptotic
or straight = HRT, DRT, JRT, ART, GRT

- ▶ For use in firm tissue, sclerotic vessels, and implants
- ▶ Needle tip triangular in section

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Blunt, round-bodied needle ○ = RN

1/2 circle, 3/8 circle or semicurved
= HRN, DRN, KRN

- ▶ For use in parenchymal tissue, the cervix, and muscle traction in the eye
- ▶ Needle tip blunt
- ▶ No puncturing of vessels or tendons

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Round-bodied needle ● = R

5/8 circle, 1/2 circle, 3/8 circle,
compound circle, J-shaped or straight
= FR, HR, DR, TR, JR, GR

- ▶ For use in soft tissue (subcutaneous) such as muscle, fascia, mucosa
- ▶ In order to improve grip by the needle holder, the cross-section of the middle part of the needle is flattened
- ▶ Easy insertion

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▶ **Absorbable sutures** approximate the tissues during the healing process. During this time the suture's tensile strength will gradually diminish. Absorbable suture material is metabolised by endogenous proteolytic enzymes or by hydrolysis (in the case of PGA *resoquick*™, GLYCOLON®, PGA RESORBA®, CAPROLON®, PDO RESORBA™ and RESORBA® Barbed Suture).

Non-absorbable sutures remain almost unchanged when placed within body tissues and are encapsulated within the wound scar tissue by the organism. The sutures used for skin

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closure are removed once the scar tissue has become sufficiently firm to hold the wound edges together (usually after 7 to 14 days).

It has to be distinguished as follows:

Absorption time

The period in which the suture loses 50 % of its knot tensile strength.

Dissolution

The period during which the suture completely dissolves within the tissue.

Disintegration

The period during which essentially non-absorbable suture breaks down by degradation into (smaller) pieces and thus losing its strength (e.g. polyamide).

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
Characteristics of absorption

Different indications also require different tensile strengths and absorption characteristics. These particular features of different sutures can be achieved by the choice of the material and modifying the production process. In addition to the immediate, moderately quick or delayed loss of tensile strength there is also the corresponding duration of absorption.

Any given thread material can only fulfil its purpose as long as it has the desired tensile strength.

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
The graph shows a simplified representation of the respective tear force losses over time. Data of GLYCOLON® and PGA RESORBA® come from in-vitro experiments, whereas the data of CAPROLON®, PDO RESORBA™ and PGA *resoquick*™ come from in-vivo experiments. Thus, the curves are not directly comparable but serve as a simplified visual representation.

RESORBA® Barbed Suture

Sterile Nadel-Faden-Kombination / needle suture combination, sterile


Nadel / needle	Stärke / gauge		Länge / length (m)	Farbe colour	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
	USP	metric				
HR, HRS = 1/2-Kreis Rundkörper HR, HRS = 1/2 circle, taper point						
HR 17 ● 	4-0	1.5	0,15	violett / violet	BSUPN0201	12
	4-0	1.5	0,30	violett / violet	BSUPN0209*	12
	3-0	2	0,15	violett / violet	BSUPN0202	12
	3-0	2	0,23	violett / violet	BSUPN0210*	12
	3-0	2	0,30	violett / violet	BSUPN0211*	12
HR 22 ● 	3-0	2	0,15	violett / violet	BSUPN0203	12
HR 27 ● 	4-0	1.5	0,15	violett / violet	BSUPN0204	12
	3-0	2	0,15	violett / violet	BSUPN0205	12
	3-0	2	0,23	violett / violet	BSUPN0206	12
	2-0	3	0,15	violett / violet	BSUPN0208	12
	2-0	3	0,30	violett / violet	BSUPN0207	12
	2-0	3	0,45	violett / violet	BSUPN0212*	12
HRS 27 ● 	2-0	3	0,15	violett / violet	BSUPN0301	12
	2-0	3	0,23	violett / violet	BSUPN0309*	12
	0	3.5	0,30	violett / violet	BSUPN0310*	12
HRS 35 ● 	2-0	3	0,30	violett / violet	BSUPN0302	12
	2-0	3	0,60	violett / violet	BSUPN0311*	12
	0	3.5	0,15	violett / violet	BSUPN0304	12
	0	3.5	0,23	violett / violet	BSUPN0305	12
	0	3.5	0,30	violett / violet	BSUPN0306	12
	1	4	0,45	violett / violet	BSUPN0303	12
HRS 40 ● 	1	4	0,45	violett / violet	BSUPN0307	12
HRS 48 ● 	1	4	0,45	violett / violet	BSUPN0308	12

RESORBA® Barbed Suture

Sterile Nadel-Faden-Kombination / *needle suture combination, sterile*

Nadel / needle	Stärke / gauge		Länge / length (m)	Farbe colour	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)	
	USP	metric					
HRN = 1/2-Kreis Rundkörper stumpf <i>HRN = 1/2 circle, blunt taper point</i>							
HRN 48 		1	4	0,45	violett / violet	BSUPN0501*	12
HRT = 1/2-Kreis schneidende Rundkörper <i>HRT = 1/2 circle, taper cutting</i>							
HRT 27 		2-0	3	0,15	violett / violet	BSUPN0601*	12
GS = Gerade schneidend <i>GS = straight, cutting</i>							
GS 50 		3-0	2	0,45	violett / violet	BSUPN1701*	12
FR = 5/8-Kreis Rundkörper <i>FR = 5/8 circle, taper point</i>							
FR 27 		3-0	2	0,23	violett / violet	BSUPN0101*	12

Healthium Medtech Limited · #472-D, 13th Cross, 4th Phase
Peenya Industrial Area, Bangalore, Karnataka-560058, India
Email: care@healthiummedtech.com
Mfg. Lic. No.: MFG/MD/2021/000367

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RESORBA Medical GmbH