PARASORB NEW RESODONT® Flex

Product features:

Resorbable in 3 – 4 months High tear strength Biocompatible Easy adaptation to the defect Does not stick to instruments




PARASORB RESODONT® Flex

PARASORB RESODONT[®] Flex is a very tear-resistant, yet flexible, resorbable barrier membrane of porcine collagen origin. Its excellent handling property makes it easy to adapt to the bony defect. PARASORB RESODONT[®] Flex does not stick to surgical instruments.

It is suitable for oral surgical procedures for guided bone and tissue regeneration.

Field of applications:

- Augmentation around dental implants
- Alveolar Ridge Augmentation
- Cover for intra-bony defects and root resection
- Periodontal Regeneration



Scanning electron micrograph of a dry PARASORB RESODONT® Flex (Cross section) at 50x magnification

REF	Size	Content
RDC1520	15 mm x 20 mm	1 membrane
RDC2030	20 mm x 30 mm	1 membrane
RDC3040	30 mm x 40 mm	1 membrane



Collagen Matrix, Inc.

15 Thornton Road, Oakland, NJ 07436, USA · **Tel:** 201.405.1477 **Fax:** 201.405.1355 · **www.collagenmatrix.com**

RESORBA Medical GmbH · Am Flachmoor 16, 90475 Nürnberg, Germany Tel. +49 9128 / 91 15 0 · Fax +49 9128 / 91 15 91 · www.resorba.com



2019-06 P1397



© 2018. Resorba is a registered trademark and company of Advanced Medical Solutions Ltd.