


RESORBA® HIGH-END COLLAGENS FOR COMMON ORAL SURGERY



YOUR DAILY
TREATMENT SUPPORT


 **RESORBA®**

REPAIR AND
REGENERATE

RESORBA is a registered trademark and company of Advanced Medical Solutions Ltd.

RESORBA® COLLAGENS

FOR COMMON ORAL SURGERY



RESORBA® is a German manufacturer and leader in the development of high-end collagens for surgical and dental wound care, as well as surgical sutures for all medical disciplines.

The generic term COMMON ORAL SURGERY describes multiple procedures that support general practitioners in their daily routine, including SOCKET PRESERVATION, an adequate TREATMENT OF PERIODONTAL INFECTIONS, MANAGEMENT OF DENTAL TRAUMAS, etc. In each case an optimal treatment is mandatory in order to reduce delayed complications above all for the patient but also for the surgeon in view of further procedures.

RESORBA® SPONGES & CONES

WITH Antibiotics

- PARASORB® Fleece Genta HD
- PARASORB® Cone Genta

WITHOUT Antibiotics

- PARASORB® Fleece
- PARASORB® Fleece HD
- PARASORB® Cone



WHY COLLAGEN?

Collagen is completely absorbable.

This means that there is no need for a second operation such as is required with non-resorbable materials. The naturally-structured collagen fibrils activate clotting like endogenous collagen.

WHY GENTAMICIN?

Gentamicin sulphate is an aminoglycoside and has a broad spectrum of antibacterial activity.

For certain antibiotics (e.g. aminoglycosides), the highest possible serum concentration determines the extent of the bactericidal action and the duration of the postantibiotic effects.

PRODUCT CHARACTERISTICS



Collagen Sponges & Cones

- Able to absorb large amounts of fluids on account of their sponge-like structure
- Fast integration into the surrounding tissue
- Rapid and safe haemostasis
- Stabilisation of the blood clot
- Easy to cut and to adapt


RESORBA® COLLAGENS

FOR COMMON ORAL SURGERY

SPONGES

ABSORBABLE WOUND COVERAGES

and inlay as a defect filler



PARASORB® Fleece Genta HD

- Highly concentrated collagen for larger bone defects
- 1 collagen gentamicin sponge of 2.5 x 2.5 x 0.5 cm contains 35 mg native equine collagen fibrils
- 25 mg gentamicin sulphate, equivalent to 13.82 – 17.88 mg gentamicin

PARASORB® Fleece

- 1 cm² contains 2.8 mg native equine collagen fibrils

PARASORB® Fleece HD

- Highly concentrated collagen for larger bone defects
- 1 cm² contains 5.6 mg native equine collagen fibrils

CONES

ABSORBABLE CONES

for treatment after tooth extractions



PARASORB® Cone Genta

- 1 collagen-gentamicin cone contains 22.4 mg native equine collagen fibrils
- 16 mg gentamicin sulphate, equivalent to 8.80 – 11.44 mg gentamicin

PARASORB® Cone

- 1 collagen cone contains 22.4 mg native equine collagen fibrils





INDICATIONS AND PRODUCT RECOMMENDATIONS

COLLAGENS WITH ANTIBIOTIC PROTECTION

PARASORB® Cone Genta

PARASORB® Fleece Genta HD

- Treatment of an extraction socket without bony defect but obvious signs of infection
- Prevention in case of „risk patients“ (diabetics, multimorbid or immunocompromised patients, etc.)
- Treatment of defective area after debridement
- Apicoectomies
- Cystectomies
- Dental Hemisections
- Root amputations
- Removal of impacted teeth in the molar area

COLLAGENS WITHOUT ANTIBIOTIC PROTECTION

PARASORB® Cone

- Treatment of an extraction socket without bony defect and without any signs of infection
- Haemostatic agent after tooth extractions and other oral surgeries

PARASORB® Fleece

PARASORB® Fleece HD

- Common Oral Mucosa Surgery
- Haemostatic agent after tooth extractions or other oral surgeries
- Wound coverage after soft tissue grafting

PLEASE NOTE



„FLEECE“-Shaped collagens are rather suggested to be used in irregular defects / „CONE“-Shaped in regular round defect morphologies.

PARASORB® Cone is a high-density collagen product.

FLEECE HD is recommended if longer absorption time is required.

RESORBA® DENTAL PRODUCTS

DENTAL COLLAGENS

Absorbable Collagen Sponges & Cones WITH Antibiotic Protection

Product	REF	Size	Content/Box
PARASORB® Fleece Genta HD	MK25	25 x 25 mm	5 fleeces
PARASORB® Cone Genta	MK11	Ø max. 12 mm height: 16 mm	5 cones single packed
	MK10	Ø max. 12 mm height: 16 mm	10 cones

Absorbable Collagen Sponges & Cones WITHOUT Antibiotic Protection

Product	REF	Size	Content/Box
PARASORB® Fleece HD	DK8001	70 x 30 mm	5 fleeces
	DK8002	18 x 36 mm	12 fleeces
PARASORB® Fleece	DK9001	70 x 30 mm	5 fleeces
	DK1836	18 x 36 mm	12 fleeces
PARASORB® Cone	DK1010	Ø max. 12 mm height: 16 mm	10 cones

Absorbable Collagen Membranes¹

Product	REF	Size	Content/Box
PARASORB RESODONT®	RD2502	22 x 25 mm	1 membrane
	RD3503	32 x 25 mm	1 membrane
	RD0703	64 x 25 mm	1 membrane
PARASORB RESODONT® Forte	RDF1502	16 x 25 mm	1 membrane
	RDF2502	22 x 25 mm	1 membrane
	RDF3503	32 x 25 mm	1 membrane
	RDF0703	64 x 25 mm	1 membrane
PARASORB RESODONT® Flex	RDC1520	15 x 20 mm	1 membrane
	RDC2030	20 x 30 mm	1 membrane
	RDC3040	30 x 40 mm	1 membrane

Manufacturer of PARASORB RESODONT® Flex:
Collagen Matrix, Inc. · 15 Thornton Road · Oakland, NJ 07436, USA

DENTAL BONE SUBSTITUTES¹

Porcine Xenograft Particulate

Product	REF	Particle Size	Content/Box
RESORBA® Xenogenic Bone Graft	XBS050	0.25 mm – 1.0 mm	1 jar, 0.5 cc
	XBS100	0.25 mm – 1.0 mm	1 jar, 1.0 cc
	XBS200	0.25 mm – 1.0 mm	1 jar, 2.0 cc
	XBS400	0.25 mm – 1.0 mm	1 jar, 4.0 cc
	XBL100	1.0 mm – 2.0 mm	1 jar, 1.0 cc
	XBL200	1.0 mm – 2.0 mm	1 jar, 2.0 cc

Manufacturer: Collagen Matrix, Inc
15 Thornton Road · Oakland, NJ 07436, USA



Biphasic Calcium Phosphate (MBCP® Technology), 60% HA / 40% β-TCP

Product	REF	Particle Size	Content/Box
RESORBA® Synthetic Bone	SB0401G50	0.5 mm – 1.0 mm	1 vial, 0.5 cc
	SB9901G01	0.5 mm – 1.0 mm	1 vial, 1.0 cc
	SB9902G02	1.0 mm – 2.0 mm	1 vial, 2.0 cc
	SB0503G05	2.0 mm – 3.0 mm	1 vial, 5.0 cc
	SB0401GS50	0.5 mm – 1.0 mm	1 syringe, 0.5 cc

Manufacturer: Biomatlante SA
ZA Les Quatre Nations · 5, Rue Edouard Belin
44360 Vigneux de Bretagne, France

Biomatlante SA is part of Advanced Medical Solutions Group plc.



¹See also our brochure for Bone Regeneration.



DENTAL SUTURES²

Absorbable Sutures

GLYCOLON® ORAL

Material Type	REF	Description	Content/Box
Polyglycolic acid-caprolactone	DB15002S	5-0 USP / 1 metric, DSM 13, 0,45 m, violet	12 pieces
Polyglycolic acid-caprolactone	DB15001S	5-0 USP / 1 metric, DSM 16, 0,45 m, violet	12 pieces
Polyglycolic acid-caprolactone	DB15003S	5-0 USP / 1 metric, DSM 18, 0,45 m, violet	12 pieces

Non-absorbable Sutures

RESOTEX® ORAL

Material Type	REF	Description	Content/Box
Polyamide	DB23001S	6-0 USP / 0.7 metric, ART 13, 0,45 m, black	12 pieces
Polyamide	DB22001S	4-0 USP / 1.5 metric, HRT 18, 0,45 m, black	12 pieces
Polyamide	DB25007S	6-0 USP / 0.7 metric, DSM 11, 0,45 m, black	12 pieces
Polyamide	DB25004S	7-0 USP / 0.5 metric, DSM 13, 0,45 m, black	12 pieces
Polyamide	DB25001S	6-0 USP / 0.7 metric, DSM 13, 0,45 m, black	12 pieces
Polyamide	DB25005S	5-0 USP / 1 metric, DSM 13, 0,45 m, black	12 pieces
Polyamide	DB25003S	5-0 USP / 1 metric, DSM 16, 0,45 m, black	12 pieces
Polyamide	DB25002S	4-0 USP / 1.5 metric, DSM 16, 0,45 m, black	12 pieces
Polyamide	DB25006S	5-0 USP / 1 metric, DSM 18, 0,45 m, black	12 pieces
Polyamide	DB28001S	4-0 USP / 1.5 metric, GR 22, 0,45 m, black	12 pieces

SUPRAMID ORAL

Material Type	REF	Description	Content/Box
Polyamide	DB32001S	4-0 USP / 1.5 metric, HRT 18, 0,45 m, black	12 pieces
Polyamide	DB31001S	4-0 USP / 1.5 metric, HS 18, 0,45 m, black	12 pieces

DENTOFIL™

Material Type	REF	Description	Content/Box
Polytetrafluoroethylene (PTFE)	DB46001	4-0 USP / 1.5 metric, HR 13, 0,45 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB42004	4-0 USP / 1.5 metric, HRT 17, 0,45 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB42003	4-0 USP / 1.5 metric, HRT 17, 0,75 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB42002	3-0 USP / 2 metric, HRT 17, 0,75 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB42001	3-0 USP / 2 metric, HRT 22, 0,75 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB45006	5-0 USP / 1 metric, DSM 11, 0,45 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB45007	5-0 USP / 1 metric, DSM 13, 0,45 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB45002	4-0 USP / 1.5 metric, DSM 13, 0,45 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB45008	5-0 USP / 1 metric, DSM 16, 0,45 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB45001	4-0 USP / 1.5 metric, DSM 16, 0,45 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB44001	3-0 USP / 2 metric, DSM 16, 0,45 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB45005	4-0 USP / 1.5 metric, DSM 19, 0,45 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB44002	3-0 USP / 2 metric, DSM 19, 0,45 m, white	12 pieces
Polytetrafluoroethylene (PTFE)	DB48001	4-0 USP / 1.5 metric, GR 22, 0,45 m, white	12 pieces

Manufacturer of DENTOFIL™:

Osteogenics Biomedical, Inc. 4620 71st Street, Bldg 78-79 · Lubbock, TX 79424 USA

²See also our brochure for Dental Sutures.



RESORBA® MORE THAN 80 YEARS OF EXPERIENCE. ARE YOU IN?



Also benefit from a world-leading independent developer and manufacturer of innovative and technologically advanced products for the global advanced wound care, surgical and wound closure markets. admedsol.com

The application pictures on the cover and on page 4 are copyrights from RESORBA/Advanced Medical Solutions Ltd. and must not be used or copied without permission.