

NYLON




- nicht-resorbierbares Material
- überdurchschnittlich geschmeidig und weich
- sehr gute Knüpfenschaften
- hervorragender Gewebedurchzug
- keine Kapillarität
- geringe Gewebereaktion

- *non-absorbable suture*
- *very soft and supple*
- *perfect handling and knotting*
- *easy passage through tissue*
- *no capillarity*
- *minimal tissue reactions*



NYLON



Eigenschaften:

- überdurchschnittlich weich und geschmeidig
- sehr gute Handhabungs- und Knüpfeigenschaften
- keine Kapillarität
- hervorragendes Gewebedurchzugsverhalten

NYLON ist ein aus Polyamid 6-6.6 extrudierter (in formbarem Zustand durch Düsen gepresster und gereckter) monofiler Faden. Aufgrund der hohen Reißfestigkeit, selbst bei geringstem Fadendurchmesser, ist Nylon speziell für feinste Nähte in der Mikrochirurgie geeignet. Polyamide können bis zu 10% Wasser binden.

Das Material ist mit einem physiologisch unbedenklichen Farbstoff eingefärbt.

- Farbe: weiß (ohne Farbstoff), schwarz
- Chemische Bezeichnung: Polyamid 6-6.6
- Fadenstärke: USP 10-0 -2 (0.2-5 metric)
- Konfektionsangebote:
 - Nadel-Faden-Kombinationen
- Sterilisation: Ethylenoxid

Properties:

- above-average softness and suppleness
- superior handling and knotting properties
- no capillarity
- excellent tissue passage

NYLON is a monofilament extruded thread (pressed and drawn through dies in a malleable condition) made from polyamide 6-6.6.







Because of its high tensile strength, even when the fibre diameter is very fine, NYLON is particularly suitable for very fine suturing in microsurgery. Polyamides can bind up to 10% water.

The material is coloured with a physiologically safe dye.

- Colour: white (no dye), black
- Chemical name: polyamide 6-6.6
- Thread diameter: USP 10-0 -2 (0.2-5 metric)
- Types of packaging:
 - needle-thread-combinations
- Sterilization method: ethylene oxide

NYLON

Sterile Nadel-Faden-Kombination / *needle suture combination, sterile*

Nadel / needle	Stärke / gauge		Länge / length (m)	Farbe / colour	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
	USP	metric				
HR F, HRS, HRX = 1/2-Kreis Rundkörper <i>HR F, HRS, HRX = 1/2 circle, taper point</i>						
2 x HR 6 F ● 	10-0	0.2	0,10	schwarz / black	5021N ¹⁾	24
HRS 35 ● 	1	4	1,50	schwarz / black	5035	24
HRS 40 ● 	0	3.5	1,50	schwarz / black	5033	24
HRX 50 ● 	0	3.5	1,50	schwarz / black	5042	24
	1	4	1,50	schwarz / black	5041	24
HRN = 1/2-Kreis Rundkörper, stumpf <i>HRN = 1/2 circle, blunt taper point</i>						
HRN 50 ○ 	1	4	1,50	schwarz / black	5050	24
HRT = 1/2-Kreis schneidende Rundkörper <i>HRT = 1/2 circle, taper cutting</i>						
HRT 18 ⊖ 	4-0	1.5	0,75	schwarz / black	5061	36

NYLON

Sterile Nadel-Faden-Kombination / *needle suture combination, sterile*

Nadel / needle	Stärke / gauge		Länge / length (m)	Farbe / colour	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
	USP	metric				
DR, DR F = 3/8-Kreis Rundkörper						
<i>DR, DR F = 3/8 circle, taper point</i>						
DR 4 ●	10-0	0.2	0,10	schwarz / black	50714N ¹⁾	24
DR 5 ●	8-0	0.4	0,15	schwarz / black	50715N ¹⁾	24
DR 5 F ●	10-0	0.2	0,15	schwarz / black	5071N ¹⁾	24
	9-0	0.3	0,15	schwarz / black	5072N ¹⁾	24
2 x DR 5 F ●	10-0	0.2	0,15	schwarz / black	5073N ¹⁾	24
DR 6 ●	10-0	0.2	0,15	schwarz / black	5074N ¹⁾	24
	9-0	0.3	0,15	schwarz / black	5075N ¹⁾	24
	8-0	0.4	0,15	schwarz / black	5076N ¹⁾	24
2 x DR 6 ●	9-0	0.3	0,15	schwarz / black	5077N ¹⁾	24
DR 10 F ●	7-0	0.5	0,45	schwarz / black	5079	36
DR 12 ●	5-0	1	0,75	schwarz / black	50712	36
DR 12 F ●	6-0	0.7	0,75	schwarz / black	50711	36
HS, HSM = 1/2-Kreis außen schneidend, mit PREMIUM-Schliff						
<i>HS, HSM = 1/2 circle, reverse cutting, PREMIUM cut</i>						
2 x HSM 8 ▽ PREMIUM-Schliff / PREMIUM cut	9-0	0.3	0,45	schwarz / black	5121N ¹⁾	24
HS 12 ▼	5-0	1	0,45	schwarz / black	5111	36
HSM 17 ▽ PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,45	schwarz / black	5126	36






















NYLON

Sterile Nadel-Faden-Kombination / *needle suture combination, sterile*

Nadel / needle	Stärke / gauge		Länge / length (m)	Farbe / colour	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
	USP	metric				
DS, DS F = 3/8-Kreis außen schneidend DSM, DSM F = mit PREMIUM-Schliff <i>DS, DS F = 3/8 circle, reverse cutting</i> <i>DS, DSM F = PREMIUM cut</i>						
DSM 6 ▽ PREMIUM-Schliff / PREMIUM cut	10-0	0.2	0,30	schwarz / black	5152N ¹⁾	24
DSM 7 ▽ PREMIUM-Schliff / PREMIUM cut	7-0	0.5	0,45	schwarz / black	5155	36
DSM 11 ▽ PREMIUM-Schliff / PREMIUM cut	7-0	0.5	0,45	schwarz / black	5156	36
	6-0	0.7	0,45	schwarz / black	5157	36
DS 13 ▼	6-0	0.7	0,45	schwarz / black	5142	36
	5-0	1	0,45	schwarz / black	5143	36
	4-0	1.5	0,45	schwarz / black	5144	36
DSM 13 ▽ PREMIUM-Schliff / PREMIUM cut	6-0	0.7	0,45	schwarz / black	51540	36
	5-0	1	0,45	schwarz / black	51513	36
	4-0	1.5	0,45	schwarz / black	51514	36
DS 16 ▼	6-0	0.7	0,45	schwarz / black	5145	36
	5-0	1	0,45	schwarz / black	5146	36
DSM 16 ▽ PREMIUM-Schliff / PREMIUM cut	4-0	1.5	0,45	schwarz / black	51517	36
DSM 16 F ▽ PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	schwarz / black	51541	36
DS 18 ▼	5-0	1	0,45	schwarz / black	5147	36
	4-0	1.5	0,45	schwarz / black	5148	36
	3-0	2	0,45	schwarz / black	5149	36
	2-0	3	0,45	schwarz / black	51423	36
	5-0	1	0,75	schwarz / black	51444	36
	4-0	1.5	0,75	schwarz / black	51445	36
	3-0	2	0,75	schwarz / black	51410	36
	4-0	1.5	0,45	weiß / white	51414	36
DS 18 F ▼	4-0	1.5	0,45	schwarz / black	51442	36
DSM 18 ▽ PREMIUM-Schliff / PREMIUM cut	5-0	1	0,45	schwarz / black	51518	36
	4-0	1.5	0,45	schwarz / black	51519	36
	3-0	2	0,45	schwarz / black	51520	36
	2-0	3	0,45	schwarz / black	51521	36
	3-0	2	0,75	schwarz / black	51522	36
DS 21 ▼	4-0	1.5	0,45	schwarz / black	51439	36
	3-0	2	0,45	schwarz / black	51440	36
DS 24 ▼	3-0	2	0,45	schwarz / black	51415	36
	3-0	2	0,75	schwarz / black	51416	36

NYLON

Sterile Nadel-Faden-Kombination / *needle suture combination, sterile*

Nadel / needle	Stärke / gauge		Länge / length (m)	Farbe / colour	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
	USP	metric				
DSM 24  ▼ PREMIUM-Schliff / PREMIUM cut	3-0	2	0,45	schwarz / black	51523	36
	2-0	3	0,45	schwarz / black	51524	36
	3-0	2	0,75	schwarz / black	51535	36
	2-0	3	0,75	schwarz / black	51536	36
DS 27 	2-0	3	0,45	schwarz / black	51419	36
DS 30 	0	3,5	0,45	schwarz / black	51421	36
	2-0	3	0,75	schwarz / black	51454	36
HSPM, HSPM F = 1/2-Kreis SPATULA, PREMIUM-Schliff <i>HSPM, HSPM F= 1/2 circle, spatulated, PREMIUM cut</i>						
HSPM 6 F  	10-0	0.2	0,15	schwarz / black	51914N ¹⁾	24
	9-0	0.3	0,15	schwarz / black	51913N ¹⁾	24
2x HSPM 6 F  	10-0	0.2	0,30	schwarz / black	5198N ¹⁾	24
	9-0	0.3	0,30	schwarz / black	5199N ¹⁾	24
2 x HSPM 7  	10-0	0.2	0,30	schwarz / black	5193N ¹⁾	24
	9-0	0.3	0,30	schwarz / black	5194N ¹⁾	24
DSPM = 3/8-Kreis SPATULA, PREMIUM-Schliff <i>DSPM= 3/8 circle, spatulated, PREMIUM cut</i>						
DSPM 4  	10-0	0.2	0,10	schwarz / black	52014N ¹⁾	24
DSPM 5  	10-0	0.2	0,30	schwarz / black	5201N ¹⁾	24
2 x DSPM 5  	10-0	0.2	0,30	schwarz / black	5202N ¹⁾	24
DSPM 6  	10-0	0.2	0,15	schwarz / black	5206N ¹⁾	24
	9-0	0.3	0,15	schwarz / black	52011N ¹⁾	24
	10-0	0.2	0,30	schwarz / black	52015N ¹⁾	24
2 x DSPM 6  	10-0	0.2	0,30	schwarz / black	5208N ¹⁾	24
2 x DSPM 7  	10-0	0.2	0,30	schwarz / black	5203N ¹⁾	24
	9-0	0.3	0,30	schwarz / black	5204N ¹⁾	24

¹⁾ Produkt hergestellt durch / product manufactured by: