



Sutures for Ophthalmological Surgery

More than 90 years of expertise
in the manufacture of surgical sutures
+ German engineering knowledge
+ Meticulous quality
= A portfolio you can rely on and trust

REPAIR & REGENERATE




Advanced Medical Solutions

RESORBA® needle-thread combinations

An unimpeded healing process is crucial for good functional and aesthetic results.

A clinician needs a highly reliable, convenient tool made of high quality materials:

PGA RESORBA®



Suggested alternative material to: VICRYL or POLYSORB

- ✓ midterm absorbable multifilament
- ✓ high tear resistance
- ✓ good handling
- ✓ excellent knot security
- ✓ atraumatic tissue passage
- ✓ minimal tissue reaction


MOPYLEN®



Suggested alternative material to: PROLENE or SURGIPRO

- ✓ non-absorbable polypropylene suture
- ✓ excellent knot security
- ✓ consistently high tear resistance
- ✓ excellent tissue passage
- ✓ hydrophobic
- ✓ non-ageing


RESOPREN®



Suggested alternative material to: PRONOVA or VASCUFIL

- ✓ non-absorbable PVDF suture
- ✓ hydrophobic, flexible and elastic
- ✓ excellent tissue passage
- ✓ extremely supple
- ✓ chemically inert
- ✓ extremely resistant to ageing

SUPOLENE



Suggested alternative material to: ETHIBOND EXCEL or TI-CRON

- ✓ non-absorbable suture
- ✓ high tear resistance
- ✓ excellent tissue passage, no sawing effect
- ✓ very even and smooth surface properties
- ✓ very minimal tissue reaction
- ✓ minimal capillarity

NYLON



Suggested alternative material to: ETHILON or MONOSOF

- ✓ non-absorbable suture
- ✓ above-average softness and suppleness
- ✓ superior handling and knotting properties
- ✓ no capillarity
- ✓ excellent tissue passage

SILK



Suggested alternative material to: PERMA-HAND SILK or SOFSILK

- ✓ non-absorbable suture
- ✓ very supple
- ✓ excellent knotting properties
- ✓ excellent knot security
- ✓ minimal sawing effect

Recommended portfolio

Selected needle-thread combinations







































	REF		Description	USP	Needle design	Thread length
▽ ◡	PAN11524 ¹⁾		PGA RESORBA® violet	7/0	DSM 7	45 cm
	7151		MOPYLEN® blue	7/0	DSM 7	45 cm
	71525		MOPYLEN® blue	6/0	DSM 7	45 cm
	4157N ¹⁾		SILK black	7/0	DSM 7F	30 cm
	4158N ¹⁾		SILK black	6/0	DSM 7F	30 cm
▽ ◡◡	4159N ¹⁾		SILK black	7/0	2 x DSM 7F	45 cm
	41510N ¹⁾		SILK black	6/0	2 x DSM 7F	45 cm
▽ ◡	PA1158		PGA RESORBA® undyed	6/0	DSM 11	45 cm
▽ ◡◡	PA11520		PGA RESORBA® violet	6/0	2 x DSM 11	45 cm
▽ ◡	771520		RESOPREN® blue	6/0	DSM 11	45 cm
	5156		NYLON black	7/0	DSM 11	45 cm
	5157		NYLON black	6/0	DSM 11	45 cm
	41517N ¹⁾		SILK black	6/0	DSM 11	45 cm
▽ ◡◡	41519N ¹⁾		SILK black	6/0	2 x DSM 11F	45 cm
▽ ◡	PA1151		PGA RESORBA® undyed	5/0	DSM 13	45 cm
	PA1152		PGA RESORBA® undyed	4/0	DSM 13	45 cm
▽ ◡◡	PA11515		PGA RESORBA® violet	5/0	2 x DSM 13	45 cm
▽ ◡	771526		RESOPREN® blue	6/0	DSM 13	45 cm
	771517		RESOPREN® blue	5/0	DSM 13	45 cm
	771527		RESOPREN® blau	4/0	DSM 13	45 cm

For a complete overview of our suture portfolio, please scan the following QR-code. You will be linked directly to our suture catalogue.



Recommended portfolio

Selected needle-thread combinations





















	REF		Description	USP	Needle design	Thread length
 	51540		NYLON black	6/0	DSM 13	45 cm
	41520N ¹⁾		SILK black	5/0	DSM 13	45 cm
	41521N ¹⁾		SILK black	4/0	DSM 13	45 cm
 	771525		RESOPREN® blue	5/0	DSM 16	45 cm
	77151		RESOPREN® blue	4/0	DSM 16	45 cm
 	51914N ¹⁾		NYLON black	10/0	HSPM 6F	15 cm
 	5198N ¹⁾		NYLON black	10/0	2 x HSPM 6F	30 cm
 	PAN1195 ¹⁾		PGA RESORBA® violet	8/0	2 x HSPM 7	30 cm
 	PA1192		PGA RESORBA® undyed	6/0	2 x HSPM 8	45 cm
	PA1193		PGA RESORBA® violet	5/0	2 x HSPM 8	45 cm
	9192		SUPOLENE weiß	5/0	2 x HSPM 8	45 cm
 	4202N ¹⁾		SILK black	8/0	DSPM 5	30 cm
 	PAN1204 ¹⁾		PGA RESORBA® violet	8/0	DSPM 6	30 cm
	PAN1205 ¹⁾		PGA RESORBA® violet	7/0	DSPM 6	45 cm
 	PAN1207 ¹⁾		PGA RESORBA® violet	8/0	2 x DSPM 6	45 cm
 	52015N ¹⁾		NYLON black	10/0	DSPM 6	30 cm
	5206N ¹⁾		NYLON black	10/0	DSPM 6	15 cm
	52011N ¹⁾		NYLON black	9/0	DSPM 6	15 cm

For a complete overview of our suture portfolio, please scan the following QR-code. You will be linked directly to our suture catalogue.



Recommended portfolio

Selected needle-thread combinations

	REF		Description	USP	Needle design	Thread length
	5208N ¹⁾		NYLON black	10/0	2 x DSPM 6	30 cm
	PA1201		PGA RESORBA® violet	7/0	2 x DSPM 7	45 cm
	PAN12014 ¹⁾		PGA RESORBA® undyed	6/0	2 x DSPM 8	25 cm
	PAN12113 ¹⁾		PGA RESORBA® violet	6/0	2 x VSPM 8	45 cm
	PA1212		PGA RESORBA® undyed	6/0	2 x VSPM 8F	45 cm
	PAN12112 ¹⁾		PGA RESORBA® violet	5/0	2 x VSPM 8	45 cm
	PA1213		PGA RESORBA® undyed	5/0	2 x VSPM 8F	45 cm
	PAN12113 ¹⁾		PGA RESORBA® violet	6/0	2 x VSPM 8	45 cm
	9211		SUPOLENE weiß	6/0	2 x VSPM 8F	45 cm
	9212		SUPOLENE weiß	5/0	2 x VSPM 8F	45 cm
	9213		SUPOLENE weiß	4/0	2 x VSPM 8F	45 cm
	4218N ¹⁾		SILK black	5/0	2 x VSPM 8F	45 cm
	7261N ¹⁾		MOPYLEN® blue	10/0	2 x GSPM 16	20 cm
	7253		MOPYLEN® blue	7/0	2 x DRTM 10	75 cm
	772510		RESOPREN® blue	6/0	2 x DRTM 12	75 cm
	772511		RESOPREN® blue	5/0	2 x DRTM 12	75 cm
	4079N ¹⁾		SILK black	4/0	DR 12	75 cm
	40615N ¹⁾		SILK black	3/0	HRT 18	45 cm

For a complete overview of our suture portfolio, please scan the following QR-code. You will be linked directly to our suture catalogue.





¹⁾ Product manufactured by: Healthium Medtech Limited · #472-D, 13th Cross, 4th Phase, Peenya Industrial Area, Bangalore, Karnataka-560058, India · Email: care@healthiummedtech.com · Mfg. Lic. No.: MFG/MD/2021/000367

RESORBA Medical GmbH

Am Flachmoor 16, 90475 Nürnberg, Germany

📍 resorba.com | ✉ infomail@resorba.com | ☎ +49 9128 / 91 15 0



Advanced Medical Solutions