

Chirurgische Federöhrnadeln steril

Spring-eyed needles sterile



Gehärteter Edelstahl:

- ✓ hohe Biegefestigkeit
- ✓ garantiert rostfrei

Optimiertes Federöhrdesign:

- ✓ sicherer Fadenhalt
- ✓ vermindert das Risiko einer Verletzung der Fadenoberfläche
- ✓ schnelles und leichtes Einfädeln durch optimale Elastizität des Federöhrs

Präzisionsschliff:

- ✓ Reduzierung des Wundtraumas durch scharfen Kantenschliff
- ✓ hohe Penetrationsfähigkeit

Tempered stainless steel:

- ✓ high bending resistance
- ✓ guaranteed rustproof










Optimised spring eye design:

- ✓ secure hold for the suture
- ✓ reduces the risk of damage to the suture's surface
- ✓ fast and easy suture insertion due to the optimal elasticity of the spring eye





Precision honing:

- ✓ sharply honed cutting edge reduces tissue trauma
- ✓ excellent penetration ability







Chirurgische Federöhrnadeln steril
Spring-eyed needles sterile

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
G (312) = 1/2-Kreis Rundkörper, standard G (312) = 1/2 circle, taper point, standard				
HR 20 ● 	G 312 /15	HR 20(70)	OESHR20N	48
HR 23 ● 	G 312 /14	HR 23(80)	OESHR23N	48
HR 26 ● 	G 312 /13	HR 26(80)	OESHR26N	48
HR 30 ● 	G 312 /12	HR 30(80)	OESHR30N	48
HR 34 ● 	G 312 /11	HR 34(90)	OESHR34N	48
HR 38 ● 	G 312 /10	HR 38(90)	OESHR38N	48
HR 42 ● 	G 312 /9	HR 42(100)	OESHR42N	48
HR 46 ● 	G 312 /8	HR 46(100)	OESHR46N	48
HR 50 ● 	G 312 /7	HR 50(110)	OESHR50N	48












Chirurgische Federöhrnadeln steril
Spring-eyed needles sterile

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
HR 54 	G 312 /6	HR 54(110)	OESHR54N	48
HR 58 	G 312 /5	HR 58(120)	OESHR58N	48
HR 63 	G 312 /4	HR 63(130)	OESHR63N	48
HR 68 	G 312 /3	HR 68(130)	OESHR68N	48
HR 73 	G 312 /2	HR 73(140)	OESHR73N	48








Chirurgische Federöhrnadeln steril
Spring-eyed needles sterile

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
E (320) = 1/2-Kreis Rundkörper, fein <i>E (320) = 1/2 circle, taper point, thin</i>				
HR 18 F ● 	E 320 /6	HR 18(60)	OESHR18FN	48
HR 20 F ● 	E 320 /5	HR 20(60)	OESHR20FN	48
HR 23 F ● 	E 320 /4	HR 23(70)	OESHR23FN	48
HR 63 F ● 	E 320 /00	HR 63(110)	OESHR63FN	48
Pb (332) = 1/2-Kreis Rundkörper, extra fein <i>Pb (332) = 1/2 circle, taper point, extra thin</i>				
HR 30 FF ● 	Pb 332 /4	HR 30(60)	OESHR30FFN	48
HR 36 FF ● 	Pb 332 /3	HR 36(60)	OESHR36FFN	48












Chirurgische Federöhrnadeln steril
Spring-eyed needles sterile

Nadel / needle		Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
GaR (314) = 1/2-Kreis Rundkörper, stark GaR (314) = 1/2 circle, taper point, strong diameter					
HR 23 S ●		GaR 314 /8	HR 23(100)	OESHR23SN	48
HR 27 S ●		GaR 314 /7	HR 27(100)	OESHR27SN	48
HR 32 S ●		GaR 314 /6	HR 32(110)	OESHR32SN	48
HR 38 S ●		GaR 314 /5	HR 38(130)	OESHR38SN	48
HR 46 S ●		GaR 314 /4	HR 46(140)	OESHR46SN	48
HR 54 S ●		GaR 314 /3	HR 54(140)	OESHR54SN	48
HR 60 S ●		GaR 314 /2	HR 60(160)	OESHR60SN	48
G (412) = 1/2-Kreis schneidende Rundkörper (Diamond Spitze) G (412) = 1/2 circle, taper cutting					
HRT 20 ▼		G 412 /15	HRT 20(70)	OESHRT20N	48
HRT 23 ▼		G 412 /14	HRT 23(80)	OESHRT23N	48
HRT 26 ▼		G 412 /13	HRT 26(80)	OESHRT26N	48
HRT 30 ▼		G 412 /12	HRT 30(80)	OESHRT30N	48









Chirurgische Federöhrnadeln steril
Spring-eyed needles sterile

Nadel / needle		Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
HRT 34 ▼		G 412 /11	HRT 34(90)	OESHRT34N	48
HRT 38 ▼		G 412 /10	HRT 38(90)	OESHRT38N	48
HRT 42 ▼		G 412 /9	HRT 42(100)	OESHRT42N	48
HRT 46 ▼		G 412 /8	HRT 46(100)	OESHRT46N	48
HRT 50 ▼		G 412 /7	HRT 50(110)	OESHRT50N	48
HRT 54 ▼		G 412 /6	HRT 54(110)	OESHRT54N	48
HRT 63 ▼		G 412 /4	HRT 63(130)	OESHRT63N	48

Chirurgische Federöhrnadeln steril
Spring-eyed needles sterile













Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
GT (212) = 1/2-Kreis außen schneidend, standard GT (212) = 1/2 circle, reverse cutting, standard				
HS 23 ▼ 	GT 212 /14	HS 23(80)	OESHS23N	48
HS 26 ▼ 	GT 212 /13	HS 26(80)	OESHS26N	48
HS 30 ▼ 	GT 212 /12	HS 30(80)	OESHS30N	48
HS 34 ▼ 	GT 212 /11	HS 34(90)	OESHS34N	48
HS 38 ▼ 	GT 212 /10	HS 38(90)	OESHS38N	48
HS 42 ▼ 	GT 212 /9	HS 42(100)	OESHS42N	48
HS 46 ▼ 	GT 212 /8	HS 46(100)	OESHS46N	48
HS 50 ▼ 	GT 212 /7	HS 50(110)	OESHS50N	48
HS 54 ▼ 	GT 212 /6	HS 54(110)	OESHS54N	48
HS 58 ▼ 	GT 212 /5	HS 58(120)	OESHS58N	48
HS 63 ▼ 	GT 212 /4	HS 63(130)	OESHS63N	48

Chirurgische Federöhrnadeln steril
Spring-eyed needles sterile

Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
HS 73 	GT 212 / 2	HS 73(140)	OESHS73N	48
GaT (214) = 1/2-Kreis außen schneidend, stark GaT (214) = 1/2 circle, reverse cutting, strong diameter				
HS 18 S 	GaT 214 / 9	HS 18(90)	OESHS18SN	48
HS 23 S 	GaT 214 / 8	HS 23(100)	OESHS23SN	48
HS 27 S 	GaT 214 / 7	HS 27(100)	OESHS27SN	48
HS 32 S 	GaT 214 / 6	HS 32(110)	OESHS32SN	48
HS 38 S 	GaT 214 / 5	HS 38(130)	OESHS38SN	48
HS 46 S 	GaT 214 / 4	HS 46(140)	OESHS46SN	48
HS 54 S 	GaT 214 / 3	HS 54(140)	OESHS54SN	48

Chirurgische Federöhrnadeln steril
Spring-eyed needles sterile



Nadel / needle	Nadelbezeichnung needle type	Euro-Code Europ. code	Best.-Nr. code no.	Packungsinhalt box content (Stück / pieces)
BT (204) = 3/8-Kreis außen schneidend, standard BT (204) = 3/8 circle, reverse cutting, standard				
DS 18 ▼ 	BT 204 /16	DS 18(70)	OESDS18N	48
DS 20 ▼ 	BT 204 /15	DS 20(70)	OESDS20N	48
DS 23 ▼ 	BT 204 /14	DS 23(80)	OESDS23N	48
DS 26 ▼ 	BT 204 /13	DS 26(80)	OESDS26N	48
DS 30 ▼ 	BT 204 /12	DS 30(80)	OESDS30N	48
DS 34 ▼ 	BT 204 /11	DS 34(90)	OESDS34N	48
DS 38 ▼ 	BT 204 /10	DS 38(90)	OESDS38N	48
DS 42 ▼ 	BT 204 /9	DS 42(100)	OESDS42N	48
DS 46 ▼ 	BT 204 /8	DS 46(100)	OESDS46N	48
DS 50 ▼ 	BT 204 /7	DS 50(110)	OESDS50N	48
DS 54 ▼ 	BT 204 /6	DS 54(110)	OESDS54N	48
DS 70 ▼ 	BT 204 /3	DS 70(130)	OESDS70N	48





RESORBA Medical GmbH

Am Flachmoor 16, 90475 Nürnberg, Germany

🖱 resorba.com | ✉ infomail@resorba.com | ☎ +49 9128 / 91 15 0



Advanced Medical Solutions